This product is disposable and is supplied sterile. Prior to the expiration date, contents are sterile if package has not been opened or damaged.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

INTENDED USE: Disposable forceps/pencils are electrosurgical devices designed to be used in soft tissue procedures. Please refer to the Kirwan Surgical Products LLC catalog or web-site (www.ksp.com) for descriptions and intended use (by surgical specialty) for individual items.

SETUP AND USE: Attach the sterile forceps/pencils to the sterile cable ensuring that the forceps/pencil pins are fully seated in the cable receptacles.

The cables to the surgical electrodes should be positioned in such a way that contact with the PATIENT or other leads is avoided. Temporarily unused ACTIVE ELECTRODES should be stored isolated from the patient.

Connect the sterile cable to the generator only while it is in the "OFF" or "STANDBY" mode. Failure to do so may result in injury or electrical shock to the patient or operating room personnel.

CAUTION: BECAUSE OF THE VARIABILITY OF OUTPUT VOLTAGES AND MODES FROM GENERATOR TO GENERATOR, DO NOT USE THE FORCEPS/PENCILS WITH GENERATORS HAVING BIPOLAR OUTPUT VOLTAGES THAT EXCEED 1100 Vp-p. REFER TO THE APPROPRIATE ELECTROSURGICAL GENERATOR MANUAL FOR INDICATIONS AND INSTRUCTIONS ON BIPOLAR OUTPUT CHARACTERISTICS TO ENSURE THAT ALL SAFETY PRECAUTIONS ARE FOLLOWED.

WARNING: CONNECT BIPOLAR ACCESSORIES TO THE BIPOLAR RECEPTACLE ONLY, AND MONOPOLAR ACCESSORIES INTO THE MONOPOLAR RECEPTACLE. IMPROPER CONNECTION OF ACCESSORIES MAY RESULT IN INADVERTENT ACCESSORY ACTIVATION OR OTHER POTENTIALLY HAZARDOUS CONDITIONS.

WARNING: DO NOT REUSE OR RE-STERILIZE. PATIENT CROSS CONTAMINATION AND/OR PRODUCT MALFUNCTION MAY OCCUR.

At the lowest power setting, test the forceps/pencils by pressing the generator's activating switch. If the generator fails to activate, check the forceps/pencils connection with the cable.

Kirwan Surgical Products LLC
180 Enterprise Drive
Marshfield, MA 02050, USA
Tel.: (781)-834-9500

Medical Product Service GmbH
Bornngasse 20
35619 Braunfels, Germany

SEE INSTRUCTIONS FOR USE

BATCH CODE

USE BY DATE

CATALOG NUMBER

SINGLE USE - DO NOT REUSE

STERILE-STERILIZED BY IRRADIATION

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