ENGLISH KIRWAN SURGICAL PRODUCTS LLC INSTRUCTIONS FOR USE AND CARE REUSABLE BIPOLAR OPHTHALMIC PENCILS

This pencil is reusable and is supplied **NON-STERILE**. Process the pencil through cleaning and sterilization prior to initial use, following guidance as outlined in this IFU. For questions or additional information on our complete line of pencils, please contact Kirwan Surgical Products LLC (KSP) at (781) 834-9500 or on the web at www.ksp.com.

INDICATIONS: The ophthalmic reusable pencil is an electrocautery device designed to be used on soft tissue surgical procedures.

WARNING: Any use of this pencil for tasks other than for which it is indicated will result in a damaged or broken pencil.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

EXTENDING PENCIL LIFE

KSP has validated these pencils for twenty-five uses. However, the number of uses obtained from the pencil depends upon the degree of care taken in processing and handling, and the surgical procedures and techniques in which the pencil is used. To extend the life of the pencil, we recommend the following:

- Not allowing gross organic contaminants to dry on the pencil (e.g., blood, mucus, and tissue), by initiating pencil decontamination immediately after completion of the surgical procedure.
- Completely drying the pencil before storage.
- Protecting the pencil from inadvertent damage while in storage, by wrapping it (especially the tip), and avoid extremes in temperatures and humidity.

NOTE: Spotting or discoloration may result from inadequate cleaning prior to sterilization or may be due to mineral deposits in water used to autoclave.

INSPECTION OF PENCIL

KSP recommends establishing a procedural review, by which the pencil is inspected frequently (before and after each use) for damage such as:

- Tip (distal end) appearance altered.
- Tip damage, e.g., burns or discoloration.
- Cracks or nicks in the molded handle of the instrument.

KSP recommends that such a review establish criteria by which a pencil showing such damage or wear is discarded and replaced.

CAUTION: Using a pencil which is damaged or worn, can be hazardous to both the patient and operating room personnel.

REPROCESSING AND STERILIZATION (i.e., cleaning & sterilization)

institutional device reprocessing and sterilization should occur in facilities that are adequately designed, equipped, monitored, and staffed by trained and experienced personnel, using the appropriate cleaning and sterilization cycle and cycle parameters. The following parameters for cleaning, and for two of the most commonly utilized methods of sterilization, are recommended as guidelines for validation.

NOTE: Reprocessing this device dictates that it undergo a thorough cleaning prior to sterilization.

WARNING: Clean and sterilize after each use.

CLEANING

- Rinse the pencil thoroughly with sterile, purified water to remove any accumulated debris.
- Hand wash the surface of the pencil (paying particular attention to the tip) using a soft bristled cleaning brush, and enzyme cleaner, (e.g., Tergazyme™ solution (Alconox, Inc.) or equivalent), to remove visible residual debris.
- **CAUTION: Avoid use of abrasive cleaners or solvents.**
- After hand washing, the surface is to be thoroughly flushed with sterile, purified water until no visible detergent residual remains.
- Once the pencil is free of cleaning solution and debris, thoroughly dry using a sterile wipe.

STERILIZATION

STEAM/GRAVITY DISPLACEMENT: DOUBLE WRAP pencil in muslin, i.e., CSR bio hospital wrap, and place (single layer) in a production tipy, steam sterilization vessel. Process at 132° C (270° F) for a 30 minute cycle.

STEAM/PRE-VACUUM: DOUBLE WRAP pencil in muslin, i.e., CSR bio hospital wrap, and place (single layer) in a production tipy, steam sterilization vessel. Process at 132° C (270° F) using pre-vacuum conditions for a 4 minute cycle.

AUTOMATED PRE-CLEANING INSTRUCTIONS

Rinse the instruments under warm running tap water until visibly clean. Use a soft bristle brush (plastic brush) as needed for hard to remove soil.

CLEANING AND DISINFECTION

Place the pencil(s) in a bath with a tested cleansing and disinfectant agent such as Renu-Klenz™ (Sliers) (14 oz/100) prepared according to manufacturer's recommendations using lukewarm tap water. The pencil(s) must be completely covered with the solution. **NOTE:** The application times, temperatures, and concentration stated by the manufacturer of the cleansing/disinfectant agent must always be observed. The pencil(s) are then immersed in the detergent solution and allowed to sonicare for ten minutes. Repeat the cleansing process if visible contamination is still present on the instrument.

Fresh solutions must be prepared daily. In case of severe soiling, the solution must be changed and replaced with a fresh solution.

A high contamination load in the ultrasonic bath impairs the cleaning action and promotes the risk of corrosion. The cleaning solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

AUTOMATED MACHINE CLEANING INSTRUCTIONS

The pencil(s) are then to be transferred via a suitable container (e.g., wire mesh basket) into the automated washer. The following cycle is recommended with these parameters programmed; set to high.

Phase	Recirculation Time (minutes)	Water Temperature	Detergent Type and Concentration
Pre-wash 1	02:00	Cold Tap water	N/A
Enzyme Wash	02:00	Hot Tap Water	Klenzyme™, 1 oz/gallon
Wash 1	02:00	65.0 C (Set Point)	Renu-Klenz™, 1/2 oz/gallon
Rinse 1	01:00	Hot Tap Water	N/A
Drying	07:00	90° C	N/A

The pencil(s) should then be dried using a clean, soft cloth and visually inspected using the naked eye under normal lighting conditions to determine that all adherent visible soil, (e.g., blood, protein substances and other debris) had been removed from all surfaces, crevices and serrations. Klenzyme and Renu-Klenz are trademarks of Sliers.

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