**ENGLISH KIRWAN SURGICAL PRODUCTS LLC**  
INSTRUCTIONS FOR USE AND CARE  
REUSABLE BIPOLAR OPHTHALMIC PENCILS

This pencil is reusable and is supplied **NON-STERILE**. Process the pencil through cleaning and sterilization prior to initial use, following guidance as outlined in this IFU. For questions or additional information on our complete line of pencils, please contact Kirwan Surgical Products LLC (KSP) at (781) 834-9500 or on the web at [www.ksp.com](http://www.ksp.com).

**INDICATIONS:** The ophthalmic reusable pencil is an electrosurgical device designed to be used in soft tissue surgical procedures.

**WARNING:** Any use of this pencil for tasks other than for which it is indicated will usually result in a damaged or broken pencil.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

**EXTENDING PENCIL LIFE**

KSP has validated these pencils for twenty-five uses. However, the number of uses obtained from the pencil depends upon the degree of care taken in processing and handling, and the surgical procedures and techniques in which the pencil is used. To achieve maximum life, KSP recommends the following:

- Not allowing gross organic contaminants to dry on the pencil (e.g., blood, mucus, and tissue), by initiating pencil decontamination immediately after completion of the surgical procedure.
- Completely drying the pencil before storage.
- Protecting the pencil from inadvertent damage while in storage, by wrapping it (especially the tip), and avoid extremes in temperatures and humidity.

**NOTE:** Spotting or discoloration may result from inadequate cleaning prior to sterilization or may be due to mineral deposits in water used to autoclave.

**INSPECTION OF PENCIL**

KSP recommends establishing a procedural review, by which the pencil is inspected frequently (before and after each use) for damage such as:

- Tip (distal end) configuration altered.
- Tip damage, e.g., burrs or discoloration.

• Cracks or nicks in the molded handle of the instrument.

KSP recommends that such a review establish criteria by which a pencil showing such damage or wear is discarded and replaced.

**CAUTION:** Using a pencil which is damaged or worn, can be hazardous to both the patient and operating room personnel.

**REPROCESSING AND STERILIZATION (i.e., cleaning & sterilization)**

Institutional device reprocessing and sterilization should occur in facilities that are adequately designed, equipped, monitored, and staffed by trained personnel. Clean and sterilize per your institution's validated procedures and cycle parameters. The following parameters for cleaning, and for two of the most commonly utilized methods of sterilization, are recommended as guidelines for validation.

**NOTE:** Reprocessing this device dictates that it undergo a thorough cleaning prior to sterilization.

**WARNING:** Clean and sterilize after each use.

**MANUAL CLEANING**

Rinse pencil thoroughly for a minimum of 30 seconds with purified/ deionized water and until visibly clean of accumulated debris.

Prepare a neutral pH enzymatic cleaner (e.g. Steris® Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner) per manufacturer's instructions.

Wash the surface of the pencil using a soft bristled cleaning brush and prepared enzymatic cleaner for a minimum of one (1) minute and until visibly clean.

**CAUTION:** Avoid use of abrasive cleaners or solvents.

Thoroughly rinse pencil by immersing it in a large volume of critical (purified) water for a minimum of 1 (one) minute. Remove the pencil and discard the rinse water. Do not reuse the water. Always use fresh volumes of water for each rinse. Repeat this stage 2 (two) additional times, for a total of 3 (three) rinses.

Once the pencil is free of cleaning solution and debris, thoroughly dry the outside surfaces with a lint free wipe and the lumen using filtered compressed air.

Visually inspect the pencil, in a well-lit area to ensure all surfaces are clean.

**STERILIZATION**

- STEAM/GRAVITY DISPLACEMENT: place pencil in a single-pouch configuration (i.e., legally approved pouch to maintain sterility) and position it (single layer) in a production type, steam sterilization vessel. Process at 132° C (270° F) for a 15 minute cycle. Dry for 20 minutes.
- STEAM/PRE-VACUUM: place forceps in a single-pouch configuration (i.e., legally approved pouch to maintain sterility) and position it (single layer) in a production type, steam sterilization vessel. Process at 132° C (270° F) using pre-vacuum conditions for a 4 minute cycle. Dry for 20 minutes.

**AUTOMATED PRE-CLEANING INSTRUCTIONS**

Rinse the instrument under warm running tap water until visibly clean. Use a soft bristle brush (plastic brush) as needed for hard to remove soil.

**CLEANING AND DISINFECTION**

Place the pencil(s) in a bath with a tested cleansing and disinfectant agent such as Renu-Klenz™ (Steris) (1/4 oz/gal) prepared according manufacturer's recommendations using lukewarm tap water. The pencil(s) must be completely covered with the solution. **NOTE:** The application times, temperatures, and concentration stated by the manufacturer of the cleansing/disinfectant agent must always be observed. The pencil(s) are then immersed in the detergent solution and allowed to sonicate for ten minutes. Repeat the cleansing process if visible contamination is still present on the instrument.

Fresh solutions must be prepared daily. In case of severe soiling, the solution must be changed sooner.

A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

**AUTOMATED MACHINE CLEANING INSTRUCTIONS**

The pencil(s) are then to be transferred via a suitable container (e.g., wire mesh basket) into the automated washer. The following cycle is recommended with these parameters programmed, set to high:

Phase	Recirculation Time (minutes)	Water Temperature	Detergent Type and Concentration
Pre-wash 1	02:00	Cold Tap water	N/A
Enzyme Wash	02:00	Hot Tap Water	Klenzyme™, 1 oz/gallon
Wash 1	02:00	65.0°C (Set Point)	Renu-Klenz™, ¼ oz/gallon
Rinse 1	01:00	Hot Tap Water	N/A
Drying	07:00	90°C	N/A

The pencil(s) should then be dried using a clean, soft cloth and visually examined using the naked eye under normal lighting condition to determine that all adherent visible soil, (e.g., blood, protein substances and other debris) had been removed from all surfaces, crevices and serrations. Klenzyme and Renu-Klenz are trademarks of Steris.

**SETUP AND USE**

**WARNING:** Connect Bipolar accessories to the Bipolar receptacle only. Improper connection of accessories may result in inadvertent accessory activation or other potentially hazardous conditions.

Connect the sterile pencil and sterile cord to the electrosurgical generator only while the generator is off (standby). Failure to do so may result in injury or electrical shock to the patient or operating room personnel.

While attaching the sterile pencil to the sterile cord ensure that the pencil pins are fully seated in the cords receptacles.

**CAUTION:** Because of the variability of output voltages and modes from generator to generator **DO NOT USE** this pencil with a Bipolar generator setting having a bipolar output voltage exceeding 1100 Vp-p. Refer to the appropriate electrosurgical generator manual for indications and instructions on voltage output characteristics to ensure that all safety precautions are followed. If no RF output is delivered to the accessory pencil when the generator's activating switch is pressed, check the cord connection with the pencil and with the generator. If proper function is still not achieved and the cord and generator function are confirmed as sound, replace the pencil and refer the questionable pencil to qualified personnel for further evaluation.

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**STERILISATION**

• **CLEAN AND STERILIZE** the pencil prior to shipping. KSP shall not accept pencil items contaminated, and a health hazard to its employees.

• Ship the pencil in a sturdy shipping box with sufficient soft packaging materials to protect it.

• Secure box with heavy tape, clearly identify the box as a return with the RGA # on the outside to expedite the process.

• Ship to:

Kirwan Surgical Products LLC  
180 Enterprise Drive  
Marshfield, Massachusetts 02050 USA  
RGA # XXXX

**Kirwan Surgical Products LLC**  For Regulatory Issues Only  
Medical Product Service  
Borgassae 20  
35619 Braunfels, GERMANY

**Manufacturer**

**REF**

**LOT**

**Date of Manufacture**

**Rx Only**

**Catalog number**

**Batch Code**

**Date of Manufacture**

**Warning, Precaution or Caution**

**Consult Instructions for Use**

**Non-Sterile**

**Not Made With Natural Rubber Latex**

**Authorized Representative in the European Community**

**KIRWAN**  
REUSABLE BIPOLAR OPHTHALMIC PENCIL  
180 Enterprise Drive, Marshfield, MA USA

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- Alıcı: Kirwan Surgical Products LLC