**ENGLISH KIRWAN SURGICAL PRODUCTS LLC  
INSTRUCTIONS FOR USE AND CARE  
REUSABLE BIPOLAR AND MONOPOLAR CORDS**

This cord is reusable and is supplied **NON-STERILE**. Process the cord through cleaning and sterilization prior to initial use, follow guidance as outlined in this IFU. For questions or additional information on our complete line of cords, please contact Kirwan Surgical Products LLC (KSP) at (781) 834-9500 or on the web at [www.ksp.com](http://www.ksp.com).

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

**INTENDED USE:** Reusable cords are designed to connect to an electrosurgical device to an electrosurgical generator. Please refer to the labeling to determine the particular connectors and specific use of the cord.

**WARNING:** Any use of this cord for tasks other than for which it is indicated will usually result in a damaged or broken cord.

**EXTENDING CORD LIFE**

KSP has validated these cords for twenty-five uses. However, the number of uses obtained from the cord depends upon the degree of care taken in processing and handling. To achieve maximum life KSP recommends the following:

- Use the molded plug at the end of the cord to disconnect. DO NOT PULL ON THE CORD'S CABLE. SEE DIAGRAM BELOW.

• Storing cords loosely coiled (4-5" diam.). Avoiding kinking, or sharply bending the cords, or placing heavy objects on them to prevent damaging the insulation or inner wire.

• Use a different cord for each procedure during the day. Keeping a supply of individually wrapped, sterile cords.

• Completely dry the cord before storing.

• Not rolling heavy tables, carts, etc., over the cord.

**INSPECTION OF CORD**

KSP recommends establishing a procedural review, by which the cord's electrical continuity is regularly tested with an ohmmeter as well as frequent inspection of the cord's insulation (before and after each use) for cracks, nicks, lacerations, or abrasions, and by which a criteria is set for the discarding and replacement of those cords which may be worn and hazardous to the patient and operating room personnel.

**WARNING:** Using a cord which is damaged, worn or until its inevitable failure, it is possible that it will overheat and either ignite itself or ignites nearby materials and is inherently dangerous to both the patient and operating room personnel.

**REPROCESSING AND STERILIZATION (i.e., cleaning & sterilization)**

Institutional device reprocessing and sterilization should occur in facilities that are adequately designed, equipped, monitored, and staffed by trained personnel. Clean and sterilize per your institutions validated procedures and cycle parameters. The following parameters for cleaning, and for five of the commonly utilized methods of sterilization, are recommended as guidelines for validation.

**NOTE:** Reprocessing this device dictates that it undergo a thorough cleaning prior to sterilization.

**WARNING:** Clean and sterilize after each use.

**MANUAL CLEANING**

Rinse Cords thoroughly for a minimum of 30 seconds with purified/ deionized water and until visibly clean of accumulated debris.

Prepare a neutral pH enzymatic cleaner [e.g. Steris® Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner] per manufacturer's instructions.

Wash the surface of the cords using a soft bristled cleaning brush and prepared enzymatic cleaner for a minimum of one (1) minute and until visibly clean.

**Caution:** Avoid use of abrasive cleaners or solvents.

Thoroughly rinse cords by immersing it in a large volume of critical (purified) water for a minimum of 1 (one) minute. Remove the cords and discard the rinse water. Do not reuse the water. Always use fresh volumes of water for each rinse. Repeat this stage 2 (two) additional times, for a total of 3 (three) rinses.

Once the cords are free of cleaning solution and debris, thoroughly dry the outside surfaces with a lint free wipe and the lumen using filtered compressed air. Visually inspect the cords, in a well-lit area to ensure all surfaces are clean.

**AUTOMATED PRE-CLEANING INSTRUCTIONS**

Rinse the instruments under warm running tap water until visibly clean. Use a soft bristle brush (plastic brush) as needed for hard to remove soil. Hard to reach areas such as, internal spaces should be flushed with a water pistol/syringe.

**CLEANING AND DISINFECTION**

Place the cord(s) in a bath with a tested cleansing and disinfectant agent such as Renu-Klenz™ (Steris) (1/4 oz/gal) prepared according manufacturer's recommendations using lukewarm tap water. The cord(s) must be completely covered with the solution. **NOTE:** The application times, temperatures, and concentration stated by the manufacturer of the cleansing/disinfectant agent must always be observed. The cords(s) are then immersed in the detergent solution and allowed to sonicate for ten minutes. Repeat the cleansing process if visible contamination is still present on the instrument.

Fresh solutions must be prepared daily. In case of severe soiling, the solution must be changed sooner.

A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visually apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

**AUTOMATED MACHINE CLEANING INSTRUCTIONS**

The cord(s) are then to be transferred via a suitable container (e.g., wire mesh basket) into the automated washer. The following cycle is recommended with these parameters programmed; set to high.

Phase	Recirculation Time (minutes)	Water Temperature	Detergent Type and Concentration
Pre-wash 1	02:00	Cold Tap water	N/A
Enzyme Wash	02:00	Hot Tap Water	Klenzyme™, 1 oz/gallon
Wash 1	02:00	65.0°C (Set Point)	Renu-Klenz™, 1/4 oz/gallon
Rinse 1	01:00	Hot Tap Water	N/A
Drying	07:00	90°C	N/A

The cords(s) should then be dried using a clean, soft cloth and visually examined using the naked eye under normal lighting condition to determine that all adherent visible soil, (e.g., blood, protein substances and other debris) had been removed from all surfaces and crevices. Klenzyme and Renu-Klenz are trademarks of Steris.

**WARNING:** Kirwan reusable cords connect electrosurgical generators to devices. Damaged cords or cords whose connectors have not been thoroughly rinsed and dried may cause electrical burns to the patient or operating room personnel.

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- STEAM/PRE-VACUUM: place cords in a single-pouch configuration (i.e., legally approved pouch to maintain sterility) and position it (single layer) in a production type, steam sterilization vessel. Process at 132°C (270°F) for a 15 minute cycle. Dry for 20 minutes.

- ETHERYLENE OXIDE (EO): DOUBLE WRAP cord in muslin i.e., CSR blue hospital wrap, and place (single layer) in a production type, EO sterilizer. Process at a nominal 600 mg/L EO concentration using Oxylume 2000 (10:90) gas for a full 2 hour cycle. Immediately following the exposure cycle aerate for 18 hours at 50°C (122°F).

- FLASH – STEAM/GRAVITY DISPLACEMENT, UNWRAPPED: Process at 134°C (273°F) for a 10 - 18 minute cycle.

- FLASH – PRE-VACUUM, UNWRAPPED Process at 132°C-134°C (270°F-273°F) for 3-18-minute cycle.

**SETUP AND USE**

Attach the sterile cord to the sterile device ensuring that the cord connector is fully seated against the device connector.

**WARNING:** Connect cords, adapters and accessories to the electrosurgical generator only while the generator is off (standby). Failure to do so may result in injury or electrical shock to the patient or healthcare provider.

**WARNING:** Connect Bipolar accessories to the **Bipolar** receptacle and **Monopolar** accessories into the **Monopolar** receptacle only. Improper connection of accessories may result in inadvertent accessory activation or other potentially hazardous conditions.

**CAUTION:** Because of the variability of output voltages and modes from generator to generator, **DO NOT use this Bipolar cord with generator setting having a Bipolar output voltage exceeding 1200Vp-p or this Monopolar cord with generator setting having a Monopolar output voltage exceeding 7000Vp-p**. Refer to the appropriate electrosurgical generator manual for indications and instructions on voltage output characteristics to ensure that all safety precautions are followed. If no RF output is delivered to the accessory handpiece when the generator's activating switch is pressed, check the cord connection with the device and with the generator. If proper function is still not achieved and the accessory handpiece and generator function are confirmed as sound, replace the cord and refer the questionable cord to qualified personnel for further evaluation.

At the lowest power setting, test the cord by pressing the generator's activating switch. **CAUTION:** The device cord should be positioned in such a way that contact with the PATIENT or other leads is avoided. Temporarily unused ACTIVE DEVICES should be stored isolated from the patient.

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**Kirwan Surgical Products LLC EC REP** Sadece Mevzuatla İlgili Sorunlar İçin