**KIRWAN SURGICAL PRODUCTS LLC**  
**INSTRUCTIONS FOR USE**  
**DISPOSABLE MONOPOLAR SUCTION COAGULATOR**

This product is disposable intended for SINGLE USE ONLY and is supplied STERILE. Prior to the expiration date, contents are sterile if the package has not been opened or damaged.

**INTENDED USE:** Suction Coagulators are electrosurgical devices designed to be used in soft tissue surgical procedures that require slow to rapid fluid evacuation and energy output for coagulation of the tissue i.e., endoscopic, sinusoidal, ENT, OB-GYN and plastic surgery.

**CONTRAINDICATIONS:** Any use of this instrument for tasks other than for which it is indicated, will usually result in a damaged or broken instrument.

**SETUP AND USE**

**WARNING:** CONNECT MONOPOLAR ACCESSORIES TO THE MONOPOLAR RECEPTACLE ONLY. IMPROPER CONNECTION OF ACCESSORIES MAY RESULT IN INADVERTENT ACCESSORY ACTIVATION OR OTHER POTENTIALLY HAZARDOUS CONDITIONS.

Attach the sterile cable to the sterile instrument ensuring that the cable connector is fully seated against the instrument connector.

**MALLEABILITY:** Shape the suction tube gradually while the stylet is inserted to help avoid its kinking.

**WARNING:** CONNECT ADAPTERS AND ACCESSORIES TO THE ELECTROSURGICAL GENERATOR ONLY WHILE THE GENERATOR IS OFF (STANDBY). FAILURE TO DO SO MAY RESULT IN INJURY OR ELECTRICAL SHOCK TO THE PATIENT OR OPERATING ROOM PERSONNEL.

**WARNING:** DO NOT REUSE OR RE-STERILIZE. PATIENT CROSS CONTAMINATION AND/OR PRODUCT MALFUNCTION MAY OCCUR.

**CAUTION:** There is a risk of gas embolism when using this device.

**CAUTION:** Activating the electrosurgical generator simultaneously with the aspiration mode may alter the path of the electrical energy away from the target.

**CAUTION:** Because of the variability of output voltages and modes from generator to generator, do not use the suction coagulators with generator settings having monopolar output voltages exceed 7000 Vp-p. Refer to the appropriate electrosurgical generator manual for indications and instructions on monopolar output characteristics to ensure that all safety precautions are followed.

Attach the suction port of the coagulator to an appropriate wall suction connector via suction tubing. Activate suction by placing gloved finger over the port in the handpiece. Keep lumen clean, use enclosed stylet. If coagulation seems diminished, scrub the tip with a soft instrument wipe in order to remove any coagulated tissue that has become adhered.

Power setting guidelines may vary due to differences in surgical techniques, patients, electrodes and surgical set-up. Start at the lowest possible power setting and increase as necessary to achieve the desired clinical effect. If the generator fails to activate, check the coagulator connection with the cable and the cable connection with the generator. If activation is still not achieved, check the switching mechanism.

 CONSULT INSTRUCTIONS FOR USE

 CATALOG NUMBER

 BATCH CODE

 USE BY DATE

 SINGLE PATIENT USE

 DO NOT RESTERILIZE

**NO LATEX** NOT MADE WITH NATURAL RUBBER LATEX

 DO NOT USE IF PACKAGE IS DAMAGED

 STERILE - STERILIZED BY IRRADIATION

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician

 Rx Only

Kirwan Surgical Products LLC  
 180 Enterprise Drive  
 Marshfield, MA 02050, USA  
 Tel: (781) 834-9500 www.ksp.com  
 Made in USA

Medical Product Service GmbH  
 Borgasse 20  
 35619 Braunfels, GERMANY

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