# KIRWAN SURGICAL PRODUCTS LLC REUSABLE COAGULATION ACCESSORIES INSTRUCTIONS FOR USE AND CARE

ENGLISH (EN).....	2
ČESKY (CS).....	5
DANSK (DA).....	8
DEUTSCH (DE).....	11
Ελληνικά (EL).....	14
ESPAÑOL (ES).....	17
SUOMI (FI).....	20
FRANÇAIS (FR).....	23
ITALIANO (IT).....	26
LIETUVIŠKAI (LT).....	29
LATVISKI (LV).....	32
NEDERLANDS (NL).....	35
NORSK (NO).....	38
POLSKI (PL).....	41
PORTUGUÊS (PT).....	44
SVENSKA (SV).....	47
TÜRKÇE (TR).....	50

## INSTRUCTIONS FOR USE AND CARE REUSABLE COAGULATION ACCESSORIES

This product is reusable, packaged in a plastic carrier (i.e. plastic tube or plastic bag) and supplied **NON-STERILE**. Process through cleaning and sterilization prior to initial use, following guidance as outlined in this IFU. For questions or additional information on our complete line of reusable coagulation accessories, please contact Kirwan Surgical Products (KSP) at (781) 834-9500 or on the web at [www.ksp.com](http://www.ksp.com).

**INTENDED USE:** Reusable coagulation handpieces are electrosurgical devices designed to be used in soft tissue surgical procedures. Reusable cords are intended to connect an electrosurgical device to an electrosurgical generator. Electrosurgical adapters are intended to convert small pins into standard banana pins or a Bovie-type pin.

**TARGET USER/PATIENT GROUP:** There is no target patient group for this device; the device is used at the surgeon's preference.

**INDICATIONS FOR USE:** Reusable coagulation handpieces are indicated for use when coagulation of soft tissue is required. Reusable cords are indicated for use when the connection of a coagulation handpiece to an electrosurgical generator is required. Electrosurgical adapters are indicated for use when the small pin attachment of a cord needs to be converted into standard banana pins or a Bovie-type pin.

**CONTRAINDICATIONS:** None known.

**CLINICAL BENEFITS:** The intended clinical benefit of this device is the effective coagulation of soft tissue.

**DEVICE CHARACTERISTICS:** These devices deliver energy from an electrosurgical generator to coagulate soft tissue.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

**CAUTION:** Any use of this device for tasks other than its intended purpose will usually result in a damaged or broken device.

**CAUTION:** Devices which are damaged or worn can be hazardous to both the patient and operating room personnel.

**CAUTION:** Do not use abrasive cleaners or solvents.

**CAUTION:** The cord to the surgical electrodes should be positioned in such a way that the contact with the PATIENT or other leads is avoided. Temporarily unused ACTIVE ELECTRODES should be stored isolated from the patient.

**CAUTION:** Because of the variability of output voltages and modes from generator to generator, **DO NOT USE** this accessory with generator setting exceeding the following voltages:

- ◆ 1100Vp-p bipolar output for bipolar pencils and bipolar adapters.
- ◆ 1200Vp-p bipolar output for bipolar forceps and bipolar cords.
- ◆ 5000Vp-p monopolar output for monopolar forceps.
- ◆ 7000Vp-p monopolar output for monopolar cords and monopolar adapters.

Refer to the appropriate electrosurgical generator manual for indications and instructions on voltage output characteristics to ensure that all safety precautions are followed. If no RF output is delivered to the accessory handpiece when the generator's activating switch is pressed, check the cord connection with the device and with the generator. If proper function is still not achieved and the accessory handpiece and generator function are confirmed as sound, replace the cord and refer the questionable cord to qualified personnel for further evaluation.

### WARNINGS AND PRECAUTIONS:

1. Clean and sterilize after each use. Use of an unsterilized device may result in infection to the patient.
2. Connect cords, adapters and accessories to the electrosurgical generator only while the generator is off (standby). Failure to do so may result in injury or electrical shock to the patient or healthcare provider. For irrigating forceps, KSP recommends gravity fed irrigation with appropriate solution. A stylet is provided for clearance of the irrigating lumen, as needed.
3. Connect **BIPOLAR** accessories to the **BIPOLAR** receptacle only, and **MONOPOLAR** accessories to the **MONOPOLAR** receptacle. Improper connection of accessories may result in inadvertent accessory activation or other potentially hazardous conditions.
4. Power setting guidelines may vary due to differences in surgical techniques, patients, electrodes and surgical set-up. Start at the lowest power setting and increase as necessary to achieve the desired clinical effect.
5. Proper disposal of devices and sharps possibly contaminated with blood, tissue, or other potentially infectious material present a biological risk and must be discarded in a closable, leak-proof, puncture-resistant receptacle, that is adequately labeled (e.g., color coding or symbology) for easy identification as biohazardous waste. Improper disposal may lead to infection of users or patients who come into contact with the device.

**EXTENDING PRODUCT LIFE:** KSP has validated these accessories for twenty-five (25) uses. However, the number of uses obtained from the accessories depends upon the degree of care taken in processing and handling, and the surgical procedures and techniques in which the accessories are used. To achieve maximum life, KSP recommends the following:

Cords	Handpieces
<ul style="list-style-type: none"> <li>• Use the molded plug at the end of the cord to disconnect.</li> <li>• Storing cords loosely coiled, avoiding sharply bending the cords, or placing heavy objects on them to prevent damaging the insulation or inner wire.</li> <li>• Use a different cord for each procedure during the day. Keeping a supply of individually wrapped sterile cords.</li> <li>• Not rolling heavy tables, carts, etc., over the cord.</li> </ul>	<ul style="list-style-type: none"> <li>• Not allowing gross or organic contaminants to dry on the accessories (e.g., blood, mucus, and tissue), by initiating accessory decontamination immediately after completion of the surgical procedure.</li> <li>• Do not allow handpiece tips to touch each other while Electrosurgical Unit (ESU) is active or activated.</li> <li>• Completely drying accessory before storage.</li> <li>• Protecting accessory from inadvertent damage while in storage, by wrapping them and avoiding extremes in temperature and humidity.</li> <li>• Avoiding STERRAD<sup>®</sup>100S on gold-plated devices.</li> </ul>

**INSPECTION OF PRODUCTS:** KSP recommends establishing a procedural review, by which the products are inspected frequently (before and after each use) for damage and to ensure products showing damage or wear are either sent for refurbishing, discarded and/or replaced. Damage to look for includes:

Cords	Handpieces
<ul style="list-style-type: none"> <li>• Cord's electrical continuity is regularly tested with an ohmmeter.</li> <li>• Inspection of the cord's insulation for cracks, nicks, lacerations or abrasions.</li> </ul>	<ul style="list-style-type: none"> <li>• Tip misalignment.</li> <li>• Tip damage e.g., burrs, bending, or discoloration.</li> <li>• For insulated instruments: cracks, nicks, lacerations, or abrasions in insulation or molded handle.</li> <li>• Cracks or nicks in the base of the instrument.</li> </ul>

**REPROCESSING** (i.e., cleaning and sterilization): Institutional device sterilization and reprocessing should occur in facilities that are adequately designed, equipped, monitored and staffed by trained personnel. Sterilize and clean per your institution's validated procedures and cycle parameters after each use. The following parameters for cleaning and the most commonly utilized methods of sterilization are recommended as guidelines for validation.

**NOTE:** Reprocessing this device dictates that it undergoes a thorough cleaning prior to sterilization.

**MANUAL CLEANING**

1. Rinse device thoroughly for a minimum of 30 seconds with purified/deionized water and until visibly clean of accumulated debris.
2. Prepare a neutral pH enzymatic cleaner (e.g. Steris® Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner) per manufacturer's instructions.
3. Wash the surface of the product using a soft bristled cleaning brush and prepared enzymatic cleaner for a minimum of one (1) minute and until visibly clean.
  - ◆ For irrigating forceps, also flush irrigating lumen with a minimum of 1mL (with prepared enzymatic detergent using a syringe.) Repeat this stage two (2) additional times, for a total of three (3) flushes.
4. Thoroughly rinse product by immersing it in a large volume of critical (purified) water for a minimum of one (1) minute. Remove the product and discard the rinse water. Do not reuse the water. Always use fresh volumes of water for each rinse. Repeat this stage two (2) additional times, for a total of three (3) rinses.
  - ◆ For irrigating forceps, also flush irrigating lumen with a minimum of 1mL with critical (purified) water using a syringe until no visible detergent remains. Always use fresh volumes of water for each flush. Repeat this stage two (2) additional times, for a total of three (3) flushes.
5. Once the product is free of cleaning solution and debris, thoroughly dry the outside surfaces with a lint free wipe and the lumen using filtered compressed air.
6. Visually inspect the product, in a well-lit area to ensure all surfaces are clean.
  - ◆ For irrigating forceps, inspect the lumen.

**AUTOMATED PRE-CLEANING INSTRUCTIONS:** Rinse the instruments under warm running tap water until visibly clean, use a soft bristle brush (plastic brush) as needed for hard to remove soil. Hard to reach areas such as, internal spaces should be flushed with a water pistol/syringe. Irrigating forceps are provided with a stylet to clean the lumen during rinsing.

**CLEANING AND DISINFECTION**

1. Place the product in a bath with a tested cleansing and disinfectant agent such as RenuKlenz™ (Steris) (1/4oz/gal) prepared according to manufacturer's recommendations using lukewarm tap water. The accessory must be completely covered with the solution.
  - ◆ The lumen of the irrigating forceps must then be flushed with the prepared detergent.
2. The product (particularly irrigating forceps) is then immersed in the detergent solution and allowed to sonicate for ten minutes. Repeat the cleansing process if visible contamination is still present on the instrument.
3. Fresh solutions must be prepared daily, in case of severe soiling, the solution must be changed sooner.
4. A high contamination load in the ultrasonic bath impairs the cleansing action and promotes the risk of corrosion. The cleansing solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling, in any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

**NOTE:** The application times, temperature and concentration stated by the manufacturer of the cleansing/disinfectant agent must always be observed.

**AUTOMATED MACHINE CLEANING INSTRUCTIONS:** The products are then to be transferred via a suitable container (e.g., wire mesh basket) into the automated washer. The following cycle is recommended with these parameters programmed; set to high.

Phase	Recirculation Time (minutes)	Water Temperature	Detergent Type and Concentration
Pre-wash 1	02:00	Cold Tap Water	N/A
Enzyme Wash	02:00	Hot Tap Water	Klenzyme™, 1oz/gallon
Wash1	02:00	65.0°C (Set Point)	Renu-Klenz™, ¼ oz/gallon
Rinse 1	01:00	Hot Tap Water	N/A
Drying	07:00	90.0°C	N/A

The device(s) should then be dried using a clean, soft cloth and visibly examined using the naked eye under normal lighting condition to determine that all adherent visible soil (e.g., blood, protein substances and other debris) had been removed from all surfaces, lumens, crevices, and serrations. Klenzyme and Renu-Klenz are trademarks of Steris.

**ENGLISH (EN)**  
**STERILIZATION**

	Forceps*	Cords/Adapters <sup>1</sup>	Ophthalmic Pencils/Probes
STEAM/GRAVITY DISPLACEMENT: place product in a single-pouch configuration (i.e., / legally approved pouch to maintain sterility) and position it (single layer) in a production type, steam sterilization vessel. Process at 132°C (270°F) for a 15 minute cycle. Dry for 20 minutes.	√	√	√
STEAM/PRE-VACUUM: place product in a single-pouch configuration (i.e., / legally approved pouch to maintain sterility) and position it (single layer) in a production type, steam sterilization vessel. Process at 132°C (270°F) using pre-vacuum conditions for a 4 minute cycle. Dry for 20 minutes.	√	√	√
CHEMICAL STERILIZATION: totally immerse the product in CIDEX® activated diaddehyde solution, (Johnson & Johnson Medical, Inc.), or equivalent. Expose the product with copious amounts of sterile water for a minimum of one minute three separate times.	√		
FLASH-STEAM/GRAVITY DISPLACEMENT, <b>UNWRAPPED</b> : Process at 134°C (273°F) for a 10-18 minute cycle.	√	√	
FLASH- PRE-VACUUM, <b>UNWRAPPED</b> : Process at 132°C-134°C (270°F-273°F) for 3-18 minute cycle.	√	√	
STERRAD® 100S: <b>DOUBLE WRAP</b> product with Spunguard® Heavy Duty Sterilization Wrap (Kimberly-Clark), or equivalent. Process a total exposure time of 50 minutes diffusion and 15 minute plasma	√		
ETHYLENE OXIDE (EO) <b>DOUBLE WRAP</b> cord in muslin i.e., CSR blue hospital wrap, and place (single layer) in production type, EO sterilizer. Process at a nominal 600mg/L EO concentration using Oxyfume 2000(10:90) gas for a full 2 hour cycle. Immediately following the exposure cycle aerate for 18 hours at 50°C(122°F).		√	

\*Remove stylet if supplied, avoid contact between cords, plug connectors, and any metal object (i.e. another instrument or tray).

<sup>1</sup>Cords whose connectors have not been thoroughly rinsed and dried may cause electrical burns to the patient or operating room personnel.

**NOTE:**

- ◆ STERRAD® 100S will discolor gold-plated forceps use another sterilization method to avoid.
  - ◆ *When sterilizing irrigating forceps, a STERRAD® 100S booster must be fitted to the irrigating lumen.*
  - ◆ Spotting/discoloration may result from inadequate cleaning or mineral deposits in water used to autoclave.
- SETUP AND USE:** Attach the sterile accessory to the sterile cord ensuring that the accessory pins are fully seated in the cord's receptacles. This condition ensures that the connection is splash proof.

**RETURNS:** For any inquiries, complaints, or serious incidents related to the use of this device, contact your Kirwan Surgical Products (KSP) distributor or KSP directly in the event a product has to be returned for evaluation. Call KSP customer service at (781) 834-9500 for a Return Merchandise Authorization (RMA). All serious incidents should also be reported to the competent authority for medical devices of your State. When returning products, follow these instructions:

- ◆ **CLEAN AND STERILIZE** prior to shipping. KSP shall not accept product it deems contaminated, and a health hazard to its employees.
- ◆ Ship product in a sturdy shipping box, with sufficient soft packaging materials to protect them.
- ◆ Secure with heavy tape, clearly identify the box as a return with the RMA# on the outside to expedite the process.
- ◆ Ship to: Kirwan Surgical Products LLC  
180 Enterprise Drive  
Marshfield, MA 02050 USA  
RMA#XXXX

Kirwan Surgical Products LLC  
180 Enterprise Drive  
Marshfield, MA 02050  
(781) 834-9500

For Regulatory Issues Only  
Medical Product Service GmbH  
Borgasse 20  
35619 Braunfels, GERMANY

Manufacturer

Catalog Number

Batch Code

Date of Manufacture

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Caution

Consult Instructions for Use

Non-Sterile

Not Made With Natural Rubber Latex

Medical Device

Authorized Representative in the European Community

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- Alıcı: Kirwan Surgical Products LLC