KIRWAN SURGICAL PRODUCTS LLC REUSABLE FORCEPS STERILIZATION TRAY INSTRUCTIONS FOR USE AND CARE

ENGLISH (EN).....	2
ČESKY (CS)	3
DANSK (DA).....	4
DEUTSCH (DE).....	5
Ελληνικά (EL)	6
ESPAÑOL (ES).....	7
SUOMI (FI).....	8
FRANÇAIS (FR).....	9
ITALIANO (IT).....	10
LIETUVIŠKAI (LT)	11
LATVISKI (LV).....	12
NEDERLANDS (NL)	13
NORSK (NO)	14
POLSKI (PL).....	15
PORTUGUÊS (PT).....	16
SVENSKA (SV).....	17
TÜRKÇE (TR).....	18

INSTRUCTIONS FOR USE AND CARE REUSABLE FORCEPS STERILIZATION TRAY

This product is reusable and packaged **NON-STERILE** within a plastic bag. For questions or additional information on our complete line of bipolar accessories, please contact Kirwan Surgical Products (KSP) at (781) 834-9500 or on the web at www.ksp.com.

INTENDED USE: The enclosed sterilization tray is designed to hold reusable bayoneted forceps during autoclave sterilization.
TARGET USER/PATIENT GROUP: There is no target patient group for this device; the device is used at the discretion of the medical center that possesses it.

INDICATIONS FOR USE: Reusable forceps sterilization trays are indicated for use when reusable bayoneted forceps have been used and need to be sterilized by autoclave sterilization.

CONTRAINDICATIONS: None known.

CLINICAL BENEFITS: The intended clinical benefit of this device is the effective support of reusable bayoneted forceps during autoclave sterilization.


DEVICE CHARACTERISTICS: This device holds reusable bayoneted forceps during autoclave sterilization.

CAUTION: Any use of this device for tasks other than its intended purpose will usually result in a damaged or broken device.
CAUTION: High pH detergents may erode anodized finish of the sterilization tray. In order to maintain the appearance of the tray's finish, a neutral pH detergent (pH 6.5-8.0) is recommended.


PRECAUTION:

1. Ensure the correct insert is inserted in the tray for the forceps intended to be sterilized. This tray comes pre-installed with the non-irrigation insert. For irrigation forceps, replace the insert in the base, where the forceps cup is received. Failure to use the correct insert may result in a damaged or broken forceps. Please refer to the images below for insert styles.

IRRIGATION

NON-IRRIGATION

SETUP AND USE:

Ensure that the tray is clean prior to use. If not, clean per standard hospital protocol.

Examine the contents of the tray itself to see that all parts are present and in working order (latches, locks, retainers, pads). If the retainers or pads have worked free from the tray, press back into position. If the locks have become loose, retighten. If any parts are missing or broken, do not use to sterilize forceps as this may cause damage to your forceps.

Disconnect the forceps to be sterilized from any connecting cables. Place forceps to be sterilized in the bottom tray. They should be adequately contained within the retainers and the tips should rest on top of the retaining pad. Do not place forceps on top of each other or allow them to contact other forceps during the sterilization process.

If the retainers will not correctly hold your device, do not use the tray for sterilization as this may cause damage to your device. Place the tray cover over the tray bottom. Depending on your tray, engage the locks or the latches at both ends of the tray. Sterilize instruments per their specific requirements or standard procedure.

For irrigation forceps, replace the insert in the base, where the forceps cup is received.

Please refer to the Kirwan Surgical Products LLC catalog for our complete line of electro-surgical devices and for replacement parts availability.

Kirwan Surgical Products LLC
180 Enterprise Drive
Marshfield, MA 02050
(781) 834-9500

Manufacturer

Date of Manufacture

Medical Device

Not Made with Natural Rubber Latex

Catalog Number

Consult Instructions for Use

Non-Sterile

For Regulatory Issues Only
Medical Product Service GmbH
Borgasse 20
35619 Braunfels, GERMANY

Batch Code

Caution

Authorized Representative in the
European Community

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