




Disposable Bipolar Suction Coagulators



Kirwan Surgical Products' line of Bipolar Suction Coagulators comes insulated and feature Kirwan's patented AURA™ Elite non-stick technology.

Kirwan's layered tube design is malleable and can withstand the rigors of electrosurgery. The gold plated tips of these coagulators allow heat to be diffused around the entire surface, eliminating hot spots and allowing for better control of the bipolar energy for more precise coagulation while maintaining a non-stick effect.

Available in 15 total styles, these coagulators are designed for soft tissue electrosurgical procedures. The 14 soft tissue coagulators feature a malleable tube, while the 14-9020 coagulators feature more rigid tubing.

ITEM #	TOTAL LENGTH	DESCRIPTION
14-9000	7 1/2" (19.1cm)	8 French
14-9001	7 1/2" (19.1cm)	11 French
14-9002	5 1/2" (14cm)	13 French
14-9004	5 1/2" (14cm)	15 French
14-9006	7 1/2" (19.1cm)	13 French
14-9008	7 1/2" (19.1cm)	15 French
14-9020*	5 1/2" (14cm)	15 1/2 French, Angled/Chiseled Tip
14-9030	11" (27.9cm)	8 French
14-9031	11" (27.9cm)	11 French
14-9032	11" (27.9cm)	13 French
14-9034	11" (27.9cm)	15 French
14-9041	7 1/2" (19.1cm)	8 French, Angled Up
14-9042	7 1/2" (19.1cm)	8 French, Angled Down
14-9043	7 1/2" (19.1cm)	11 French, Angled Up
14-9044	7 1/2" (19.1cm)	11 French, Angled Down

*more rigid tubing