



Reusable Non-stick Forceps for Microsurgical and Plastics indications

Kirwan's reusable microsurgical forceps provide the best opportunity for excellent results with all plastic and reconstructive surgeries. The elite non-stick qualities minimize potential damage to delicate vessels at risk.

The elite non-stick qualities of these reusable forceps ensure the best possible outcome for plastics cases including: DIEP plastic reconstruction, breast reconstruction, lateral flaps, abdominoplasty, free flaps, carpal tunnel, micro hand surgery, thyroid resections, ocular plastics, and facial plastics. AURA non-stick forceps employ highly conductive metals throughout the entire body. Therefore, the tips are not coated and the non-stick effect will not diminish over time.

10-3402AU Reusable 4" Jewelers Non-stick Bipolar Forceps, Insulated with 0.4mm tips. The four inch standard Jewelers forceps has been widely used in microsurgical procedures for decades. In addition to its exceptional non-stick features the 10-3402AU's fine tips offer precise tissue coagulation where it's needed most.

10-3442AU Reusable 5" Jewelers Non-stick Bipolar Forceps, Insulated, with 0.4mm tips. A Jewelers with an elongated neck and gripped handle. This forceps has been extremely popular and allows the surgeon precise tissue coagulation in slightly deeper surgical sites while not sacrificing comfort and performance.

10-3415AU Reusable 5" Adson Non-stick Bipolar Forceps, Insulated, with 1.0mm tips. The Adson's comfortable grip and broad tip allow for greater heat dispersion and an excellent non-stick effect even when used in fatty tissue.



180 Enterprise Drive Marshfield, MA 02050
Phone: 781-834-9500 • 888-547-9267 • Fax: 781-834-0022
Web: www.ksp.com • Email: info@ksp.com